

MTXQCvX2 documentation

Christin Zasada

2019-02-07

Contents

| | | |
|-----------|--|-----------|
| 1 | Welcome | 5 |
| 2 | Content of MTXQCvX documentation | 7 |
| 3 | Workflow for Maui-annotation projects | 9 |
| 3.1 | Read this in case | 9 |
| 3.2 | Quick view | 9 |
| 3.3 | Input files | 9 |
| 3.4 | Annotation-file | 10 |
| 3.5 | Sample_extracts-file | 10 |
| 3.6 | Internal Standard | 11 |
| 4 | Workflow for Metmax-extracted projects | 13 |
| 4.1 | You want to follow this | 13 |
| 4.2 | Introduction | 13 |
| 4.3 | Quick view | 13 |
| 4.4 | Input files | 14 |
| 4.5 | Annotation-file | 14 |
| 4.6 | Sample_extracts-file | 14 |
| 4.7 | Update metabolite names in conversion_metabolite.csv | 15 |
| 5 | MTXQCvX2_init | 17 |
| 6 | MTXQCvX_experimentalSetup.Rmd | 19 |
| 7 | MTXQCvX_part1.Rmd | 21 |
| 8 | MTXQCvX_part2.Rmd | 23 |
| 9 | MTXQCvX_part3.Rmd | 25 |
| 10 | MTXQCvX_part4.Rmd - Metmax parser | 27 |
| 10.1 | This section explains | 27 |
| 10.2 | Input files | 27 |
| 11 | Configuration of MTXQCvX2 - config_mtx/... | 29 |
| 11.1 | conversion_metabolite.csv | 29 |
| 11.2 | letter_pathway_complete.csv | 29 |
| 11.3 | quant1-values.csv | 30 |
| 11.4 | incorporation_calc.csv & mid_backups.csv | 30 |
| 12 | Protocols - Sample extraction | 31 |
| 12.1 | Cell extracts | 31 |

| | |
|---|-----------|
| 12.2 Tissue samples | 32 |
| 12.3 Blood samples | 32 |
| 13 Protocols - GC-ToF-MS measurement | 33 |
| 13.1 Sample derivatisation | 33 |
| 13.2 GC-MS measurement | 33 |
| 14 Data processing - MAUI | 35 |
| 14.1 Processing In ChromaToF | 35 |
| 14.2 Maui notes | 36 |
| 14.3 Maui exports | 36 |
| 15 Data Processing - Metmax | 37 |
| 15.1 Resampling | 37 |
| 15.2 1D-basic | 37 |
| 15.3 Reference search | 37 |
| 15.4 Export for Metmax | 37 |
| 15.5 Data extraction with Metmax | 37 |
| 16 Frequently Asked Questions | 39 |
| 16.1 What are additional quantification standards | 39 |
| 16.2 How do I extend conversion_metabolite.csv | 39 |
| 17 Appendix - Tables | 41 |
| 17.1 MTXQC base tables | 41 |
| 17.2 Input data | 44 |
| 17.3 Output data | 44 |
| 18 pSIRM experiments | 51 |
| 18.1 Experimental design | 51 |
| 18.2 Experimental procedures | 52 |
| 18.3 Protocol pSIRM | 52 |
| 18.4 Hints & notes | 54 |
| 19 GC-MS settings | 55 |
| 20 Solutions and Materials | 59 |
| 20.1 Solutions | 59 |
| 20.2 Materials | 60 |

Chapter 1

Welcome

This documentation introduced to you how to use MTXQCvX2 in order to assess the quality of your GC-MS derived data, perform the determination of calibration curves and absolute quantification. It furthermore provides you two normalisation strategies and the calculation of quantities in, e.g., pmol/1e+6 cells or pmol/mg tissue.

MTXQCvX2 does also enable the calculation of stable isotopic incorporation and the evaluation of the underlying data, the mass isotopomer distributions (MIDs).

The tool has been set up to support the in-lab developed workflow for quantitative metabolomics experiments using the in-house developed software Maui for the annotation of data. MTXQCvX2 bridges the gap between quality control and first data post-processing / analysis of GC-MS derived data (MTXQCvX2_part1, MTXQCvX2_part2).

Nevertheless MTXQCvX2 includes a module in order to integrate all kind of data provided in spreadsheet-format, e.g., derived from metmax, extracting required information and creating corresponding files (MTXQCvX2_part4).

Both workflows are introduced in the distinct chapters including their required input parameter (chapter Workflow Maui and Workflow Metmax).

Chapter 2

Content of MTXQCvX documentation

This documentation introduces to you how to use MTXQCvX2 in order to run a first straight-forward data analysis of your metabolomics experiment. Experimental and mathematical concepts have been introduced for the pulsed stable isotope resolved metabolomics (pSIRM) approach published in (Pietzke et al., 2014).

MTXQCvX2 contains a suite of modules is optimized to process GC-MS derived data and processed either in Maui or Chromatof/Metmax. Workflows for both approaches are introduced with step-by-step instructions in chapter @??maui) and @??metmax).

Subsequently each MTXQCvX2 module (chapter 6 - ??) is introduced in greater detail including a input / output files and processing parameter.

The configuration of MTXQCvX2 has been split into two categories - (1) a general configuration `config_files` and (2) metabolomics specific parameters `config_mtx`. Latter one is meant to provide flexibility including further substances. How to do so and what files can be customizes is shown in chapter ??.

Protocols of sample preparation and measurement are summarised in the following chapters ?? and ??. Workflow-specific processing methods applied Chromatof are introduced separately in the chapter 14 and ?? including all parameter.

The remaining chapter cover a chapter for frequently asked questions (chapter16) and the visualisation of all tables (chapter 17).

Chapter 3

Workflow for Maui-annotation projects

3.1 Read this in case

- you have run a Maui project
- exported all required container (see ??)
- you have a copy of sequence list and experimental conditions
- you know the extraction procedure

The following article describes briefly how to use MTXQCvX2 in case you used Maui for the annotation of your metabolomics project. It does not matter if you have performed an experiment including stable isotopes or if you just aim for the quantification of a few intermediates.

3.2 Quick view

1. Setup a new R-project and copy MTXQC template files and folders
2. Knit with parameter: `MTXQC_init.Rmd` and create project folder, e.g., `psirm_glucose`
3. Copy input files and rename `ManualQuantTable.tsv` (`e18205cz.tsv`)
4. Create `annotation.csv` and `sample_extracts.csv` files
5. Define the internal extraction standard
6. Knit with parameter: `MTXQC_ExperimentalSetup.Rmd`
7. Knit with parameter: `MTXQC_part1.Rmd`
8. Knit with parameter: `MTXQC_part2.Rmd`
9. If required, proceed with `MTXQC_part3.Rmd` for ManualValidation

3.3 Input files

Three different kind of export functions have been implemented in Maui. These functions provide the export of the actual data into `.csv` or `.tsv` files that are directly usable as input files for MTXQCvX2. Please refer to section 14.3 how you perform the export and which containers have to be exported using what export function and where to copy them in `psirm_glucose/input/`.

Certain circumstances might wish you to combine *multiple MAUI-projects* into one MTXQC-project. This might be the case when you run the same samples in split and splitless mode on the machine or your

experimental setup has been measured in multiple batches in order to avoid derivatisation effects.

It is highly recommended to combine the input files derived from a different Maui projects beforehand the analysis. In that way you have only to work with a single file `CalculationFileData.csv`¹ containing all data of the your experiment.

The herein described process provides a quick way how to combine the exported files from different Maui projects. The script `combine-sets.R`² automatically saves all combined files into the correct `input` folder. Just update the folder and subfolder names. All the rest has been taken care of for you.

1. Create in the MTXQC-project folder (e.g., `psirm_glucose/`) a new folder called `raw-data`
2. Create a subfolder for each Maui-project in `psirm_glucose/raw_data/...`
3. Copy into this folder all your Maui-derived input files altogether
4. Update the parameter of `combine-sets.R`, meaning folder name definitions, file
5. Execute the R script
6. Merged files have been generated and copied into the corresponding folder: `psirm_glucose/input-folder/gc/...` or `psirm_glucose/input-folder/inc/...`
7. Copy the renamed `ManualQuantTable.tsv` files of each Maui project into `psirm_glucose/input/quant/...`

3.4 Annotation-file

The annotation file relate file names with experimental conditions or specify quantification standards in your batch. Two columns - **File and Type** - are obligatory and have to be present in the annotation file. In the case of their absence `MTXQCvX_part1.Rmd` stops processing and shows an error message.

A quick way to generate an annotation file is described below:

1. Copy the first row / header of `quantMassAreaMatrix.csv` file
2. Paste & transpose the content into a new Excel-File into column A
3. Change the first entry: Metabolite -> File
4. Remove the entry `QuantMasses` at the very end of the column A
5. Add the column `Type` and specify each file either as **sample** or **addQ1_dilution**³
6. Add more columns specifying your experimental conditions, e.g., Cellline and Treatment⁴
7. Save the content as csv-file in the `psirm_glucose/input/...`

3.5 Sample_extracts-file

The `sample_extracts.csv` file is required in order to determine automatically absolute quantities in the manner of pmol/1e+6 cells or pmol/mg tissue in the `CalculationFileData.csv`.

This file requires two obligatory columns - **Extract_vol** and **Unit**⁵. Please specify for each experimental condition the amount of extracted cells (count), tissue (mg) or volume of blood/plasme (ul) in the unit shown in the brackets.

The names of the columns of the experimental conditions have to match up with the annotation file. Save the file in the folder `psirm_glucose/input/....`

If the defined experimental conditions do not match up with the annotation `MTXQCvX2_part1.Rmd` exit data processing. A template file is saved for review and usage at `inst/template_files/...`

¹stored in `psirm_glucose/output/quant/...`

²`inst/template_files/...`

³see for further details `additionalQuant`

⁴optimal: two-three parameter, max: four parameter. Consider possible combinations, e.g., HCT116-control, HCT116-BPTES

⁵Define: count, mg or ul

3.6 Internal Standard

MTXQCvX2 allows the specification of project-specific internal extraction standards. The only thing you need to do is to define the corresponding compounds as an internal standard in the `conversion_metabolite.csv` file. To do so, add `InternalStandard` to the compound in last column `Standard`.

For an classical pSIRM experiment in the Kempa lab we are using cinnamic acid. The evaluation of this compound has been integrated in Maui. Peak areas of cinnamic acid are exported from a distinct container called `cinAcid`. The exported file has to be renamed to `InternalStandard.csv` though and moved to `psirm_glucose/input/gc/....`.

If you have used a different compound as an internal extraction standard you might need to extract the peak areas of this compound from the Maui export file `quantPeakAreasMatrix.csv` file and save it in the folder `psirm_glucose/input/gc/InternalStandard.csv`, respectively. Prerequisite - you have annotated the compound in Maui.

The report of `MTXQCvX2_part1.Rmd` includes the defined internal standard for each project in a message.

Chapter 4

Workflow for Metmax-extracted projects

4.1 You want to follow this ...

- in case you have measured samples and quantification standards by GC-MS
- performed the annotation of intermediates in ChromaToF or vendor software
- exported all information into .txt files
- used metmax to extract peak areas / mass isotopomer distributions (MIDs)

4.2 Introduction

This document describes how to use MTXQCvX2 in combination with metmax¹.

Historically, MTXQCvX2 has been developed and optimized for Maui-derived input files. The MTXQCvX2-part4.Rmd functions as a converter of metmax-derived files in order to create suitable input formats for MTXQCvX-part1.Rmd.

This module could also be used to convert tables derived from other programs as long as they are stick with the herein described table formats. Mandatory columns are referenced in the text for each kind of input file.

The general workflow of the NMTXQCvX2 project is briefly shown below in quick view. More detailed instructions are summarised in the following paragraphs.

For more detailed explanations about the individual input parameter for each module of MTXQCvX2 please proceed to read the documentation about the individual modules and their knitting parameter. The relation of knitting parameter, input and output files are described in each section.

4.3 Quick view

1. Generate input files: run MTXQC_part4.Rmd²
2. Setup R-project and copy MTXQC-files
3. Knit with parameter: MTXQC_init.Rmd
4. Copy input files into corresponding folders

¹<http://gmd.mpimp-golm.mpg.de/apps/metmax/default.htm>

²read here the instructions

5. Create annotation.csv and sample_extracts.csv files³
6. Update metabolite names in `conversion_metabolite.csv`⁴
7. Define the internal standard and/or alkanes⁵
8. Knit with parameter: `MTXQC_ExperimentalSetup.Rmd`
9. Knit with parameter: `MTXQC_part1.Rmd`
10. Knit with parameter: `MTXQC_part2.Rmd`
11. If required - proceed with `MTXQC_part3.Rmd` for ManualValidation

4.4 Input files

If you need an introduction about how to use metmax - have a look at the separate documentation `Metmax_intro`.

The chapter ?? `MTXQCvX_part4` explains in detail how to use this module to generate suitable input files.

4.5 Annotation-file

The annotation file relate file names with experimental conditions or specify quantification standards in your batch. Two columns - **File and Type** - are obligatory and have to be present in the annotation file. In the case of their absence `MTXQCvX_part1.Rmd` stops processing and shows an error message.

A quick way to generate an annotation file is described below:

1. Copy all file names from a file of your choice
2. Paste & transpose the content into a new Excel-File into column A
3. Call column A -> File
4. Optional: Remove any non-file name entry in this column
5. Add the column Type and specify each file either as **sample**, **Q1_dilution**, **addQ1_dilution**⁶
6. Add more columns specifying your experimental conditions, e.g., Cellline and Treatment⁷
7. Save the content as csv-file in the `psirm_glucose/input/...`

4.6 Sample_extracts-file

The `sample_extracts.csv` file is required in order to determine automatically absolute quantities in the manner of pmol/1e+6 cells or pmol/mg tissue in the `CalculationFileData.csv`.

This file requires two obligatory columns - **Extract_vol** and **Unit**⁸. Please specify for each experimental condition the amount of extracted cells (count), tissue (mg) or volume of blood/plasme (ul) in the unit shown in the brackets.

The names of the columns of the experimental conditions have to match up with the annotation file. Save the file in the folder `psirm_glucose/input/....`

If the defined experimental conditions do not match up with the annotation `MTXQCvX2_part1.Rmd` exit data processing. A template file is saved for review and usage at `inst/template_files/...`

³Details further down this document

⁴Column: `Metabolite_manual`

⁵Also in `conversion_metabolite.csv`; see below paragraph Standards

⁶see for further details `additionalQuant`

⁷optimal: two-three parameter, max: four parameter. Consider possible combinations, e.g., HCT116-control, HCT116-BPTES

⁸Define: count, mg or ul

4.7 Update metabolite names in `conversion_metabolite.csv`

The file `conversion_metabolite.csv`, saved in `config_mtx/`, serves as a kind of translational table. It defines alternative version of metabolite library names that come in handy to plot data using shorter metabolite names. This file is also used to define settings and standard classifications. Detailed information for each column of the file are shown here: REF

4.7.1 Match your annotation with library names

Prior the analysis you need to match the names of your intermediates with the `conversion_metabolite.csv` file. You need to update or add the corresponding name for each intermediate in the column **Metabolite_manual**.

General suggestion for naming conventions in ChromaToF: Metabolite_Derivate, e.g., Lactic acid_(2TMS). In case of the presence of main- (MP) and byproducts (BP) use: Metabolite_Derivate_MP/BP, e.g., Glucose_(1MEOX)(5TMS)_MP.

If you have annotated intermediates that are not included so far in this table please follow the instructions how to extend `conversion_metabolite.csv`.REF

4.7.2 Define your internal standards and alkanes

MTXQCvX2 allows the specification of project-specific internal standards. Corresponding compounds have to be marked as an internal standard in `conversion_metabolite.csv` by adding the tag **InternalStandard** in the column Standard.

If you check the box - InternalStandard in the parameter selection for `MTXQCvX2_part4.Rmd` the module searches in your input file for peak areas of the defined standard and extracts the information. It also generates the file `InternalStandard.csv` and stores it at `psirm_glucose/input/gc/...`

In the same way alkanes are defined in `conversion_metabolite.csv`. Each alkane has to be flag tagged with **Alk** in the column Standard. This gives you the opportunity to implement customized mixtures of alkanes in order to determine the retention index. `MTXQCvX_part4.Rmd` recognises the flag tag and generates `Alcane_intensities.csv` based on your input file containing peak areas and saves it in `psirm_glucose/input/gc/...`⁹.

The in-lab protocol considers nine alkanes from c10 to c36. Standard annotation includes an hashtag, e.g., #c10. If you use this annotation even Metmax would be able to determine the retention index.

⁹It should be alkane, I know, but Maui doesn't, unfortunately...

Chapter 5

MTXQCvX2__init

MTXQCvX2__init.Rmd - why and how to use it. Advantages of the project folder.

Chapter 6

MTXQCvX_experimentalSetup.Rmd

Chapter 7

MTXQCvX_part1.Rmd

Chapter 8

MTXQCvX_part2.Rmd

Chapter 9

MTXQCvX_part3.Rmd

Chapter 10

MTXQCvX_part4.Rmd - Metmax parser

10.1 This section explains ...

- what MTXQCvX_part4.Rmd does
- how do input files need to look like
- which files are generated
- what the distinct checkboxes mean

This module provides the generation of suitable input files for MTXQCvX2 based on spreadsheet exported information by tools like metmax.

10.2 Input files

10.2.1 Quantification - PeakAreas.csv

In order to perform absolute quantification of

You need a file containing all extracted peak areas for each metabolite and file¹. The header of metmax-extracted files looks like shown below (see table 1). Please, remember to delete the second header row, representing the column loads for each file before saving as csv-file. Otherwise you end up with weird imported dataframes in R. Quantification masses have to be updated while processing in ChromaToF prior the export of the data e.g., with a reference search² or using statistical compare. pSIRM experiments require the definition of pTop5 masses³ instead of top5 masses in the reference in order to take into account the shift of intensities induced by the application of stable isotopes⁴

| name | mass | ri | row.load | file_1 | file_2 | file_x |
|------|------|------|----------|--------|--------|---------|
| Lac | 219 | 1051 | 0.76 | 15423 | 135444 | 465486 |
| Pyr | 174 | 1042 | 0.65 | 56978 | 46888 | 4354544 |
| Cit | 273 | 1805 | 0.99 | 1326 | 23321 | 132121 |

¹Tools/Options/Retention analysis, Parameter: Area

²See `vignette/ReferenceSearch`

³Extended list of quant masses considering isotope incorporation

⁴Mandatory columns: name, mass, files

MTXQCvX_part4 takes care of the formatting and correct column names of the peak areas file and saves it⁵. MTXQCvX_part4 generates also the file PeakDensities-Chroma.csv⁶, in case you have selected the option to include sum of area normalisation while knitting this module.

10.2.2 Isotope incorporation - MIDs.csv⁷

In order to determine the incorporation of stable isotopes MTXQCvX2 requires as an input the mass isotopomer distributions (MIDs) for each intermediate and measurement⁸. Fragments for each intermediate have to be pre-defined in metmax at Tools/Options/metabolite masses. They can be imported⁹ or manually specified each by each. An example of the metmax output is shown in table 2. The output has to be saved as csv-file, including the deletion of the partial row `column.load`, respectively¹⁰.

| name | mass | ri | row.load | file_1 | file_2 | file_x |
|------|------|------|----------|--------|---------|--------|
| Lac | 219 | 1051 | 0.85 | 31026 | 5165829 | 5829 |
| Lac | 220 | 1051 | 0.85 | 3607 | 662277 | 277 |
| Lac | 221 | 1051 | 0.85 | 1222 | 111481 | 81 |
| Lac | 222 | 1051 | 0.85 | 188 | 1003494 | 10023 |
| Lac | 223 | 1051 | 0.85 | 0 | 33542 | 342 |

MTXQCvX_part4 calculates the stable isotope incorporation and exports DataMatrix.csv as well as pSIRM_SpectraData.csv¹¹. The mathematics behind are outlined in (Pietzke et al., 2014)

Important: Extracted MIDs have to match with defined mass couples for each metabolite in MTXQCvX2¹². Please refer for more details to `vignettes/config_mtx-files`.

10.2.3 Derivatisation efficiency - mz73.csv¹³

The extraction of intensities for the ion m/z 73 works analogous to the extraction of MIDs¹⁴. Mass ranges have to be defined for each intermediate for the mass 73 by defining starting and end mass with 73. MTXQCvX_part4 generates the file MassSum-73.csv¹⁵. Check `inst\template_files\` for reference. Hopefully soon a new metmax button extracting specific intensities across the batch.

⁵input/quant/quantMassAreasMatrix.csv

⁶input/gc/PeakDensities-Chroma.csv

⁷Required for calculation isotope incorporation

⁸Tools/Options/Isotope concentrator; Parameter: IntensityOfMass

⁹inst/template_files/MetMax_MIDs.txt

¹⁰Mandatory columns: name, mass, files

¹¹input/inc/DataMatrix & pSIRM_SpectraData.csv

¹²config_mtx/incorpo_calc_masses.csv

¹³Required for: sum of area normalisation

¹⁴Tools/Options/Isotope concentrator; Parameter: IntensityOfMass

¹⁵input/gc/MassSum-73.csv

Chapter 11

Configuration of MTXQCvX2 - config_mtx/...

Herein explained are the customizable tables of the MTXQCvX2 universe.

11.1 conversion_metabolite.csv

| Column.name | Description | Value |
|-------------------|---|-------------------------------------|
| Metabolite_manual | Manual defined metabolite name | #Alanine (2TMS) |
| Metabolite | Library name of the metabolite | Alanine_(2TMS)_BP_RI:1097_IDENT:B+C |
| Metabolite_short | Short version of library name of the metabolite | Alanine_(2TMS) |
| Lettercode | Lettercode version of metabolite name | Ala_2TMS |
| Q1_value | Checked if quant1:1 value available | x |
| Mass_Pos | m/z-value corresponding to m_inc | 118 |
| SE_sel | Evaluation of the MIDs | x |
| Q_sel | Evaluation for absolute quantification | x |
| nopsirm | Exclusively for absolute quantification | |
| Standards | Defined as standard | InternalStandard, Alk |

11.2 letter_pathway_complete.csv

| Column.name | Description | Value |
|-----------------|---|---------------|
| Letter_Derivate | Derivate definition | Ala |
| Lettercode | Lettercode name of metabolite | Ala_3TMS |
| Pathway | Ass.pathway | aa |
| Pathway.1 | Ass. pathway - ordered for heatmap | 5-aa |
| Met_pathway | Ass. pathway - ordered for heatmap incl. Lettercode | 5-aa_Ala_3TMS |
| Subs_class | Substance class | aa |
| Met_class | Substance class incl. Lettercode | aa_Ala_3TMS |

11.3 quant1-values.csv

| Column.name | Description | Value |
|-----------------|-----------------------------|-------|
| Letter_Derivate | Derivate name of metabolite | 3PGA |
| Quant1_v4 | Quantity in (pmol) | 43480 |
| Quant1_v3 | Quantity in (pmol) | 43480 |

11.4 incorporation_calc.csv & mid_backups.csv

| Column.name | Description | Value |
|-------------|----------------------------|-------------------------------------|
| Metabolite | Library name of metabolite | Alanine_(2TMS)_BP_RI:1097_IDENT:B+C |
| Mass_mz | m/z-value | 116, 118 |
| LI_MID | Definition of mass level | m0, minc |

| Column.name | Description | Value |
|----------------|---------------------------------------|--|
| Metabolite | Library name of metabolite | Alanine_beta_(3TMS)_MP_RI:1435_IDENT:A+D |
| Mass.m.z. | m/z value | 188 |
| BackupPeakArea | Peak area of Backup MID | 4960 |
| BackupMID | MID value for corresponding Mass.m.z. | 0.8005 |

Chapter 12

Protocols - Sample extraction

12.1 Cell extracts

Materials:

- cell culture dishes (10 cm), max. confluency 75%
- washing buffer (Hepes, NaCl, pH 7.4)
- 50% MeOH, ice-cold
- 2 mg/ml cinnamic acid
- chloroform
- 15 ml falcon tubes
- cell lifter

Procedure:

- prepare cell culture dishes accordingly to your experimental conditions
- discard cell culture media
- add quickly 5 ml of washing buffer, discard it
- add very immediately 5 ml ice-cold 50% MeOH suppl. 2 µg/ul cinnamic acid
- detach cells using cell lifter
- collect and transfer cell extract into 15 ml falcon
- store falcons until further processing on ice
- add 1 ml chloroform
- incubate for 60 min at cold temperatures (4 °C) on rotary or thermo shaker
- centrifuge at max speed for 10 min, cold temperatures
- collect polar and lipid phases into fresh falcons / tubes
- dry under vacuum

In order to generate technical backups:

- resuspend dried extracts in 600 µl 20% MeOH
- shake at cold temperature on thermo shaker for 30 min
- split volumes into equal parts in fresh eppendorf tubes
- dry under vacuum

Suggested cell density: $2 - 3 \times 10^6$ cells / extract.

12.2 Tissue samples

Materials:

- Methanol:Chloroform:Water (MCW) in ratio 5:2:1
- 2 mg/mg cinnamic acid in MeOH
- ddH₂O
- eppendorf tubes
- tissue lyzer / pulverizer

Procedure:

- snap-freeze tissue samples
- pulverize samples
- aliquote 50 mg of tissue powder
- add 1.5 ml of MCW (suppl. with cinnamic acid final conc. 2 ug/ul)
- shake for 60 min on rotary shaker at cold temperature (4 C)
- add 0.5 ml ddH₂O for phase separation
- centrifuge maximum speed, 10 min, cold temperatures
- collect polar and lipid phases in fresh vessels
- dry under vacuum

12.3 Blood samples

Material:

- Methanol:Chloroform:Water (MCW) in ratio 5:2:1
- 2 mg/mg cinnamic acid in MeOH
- ddH₂O
- eppendorf tubes

Procedure:

- give 20 ul blood / sera directly into 1 ml MCW to avoid lumps
- in case of lumps sonicate samples
- shake samples at 4 C for 800 rpm for 60 min
- add 500 ul ddH₂O and vortex shortly
- spin down at 4 C at max speed for 10 min
- aliquote polar phase into 2-3 times 500 ul in 1.5 ml tubes
- aliquote lipid phase 2 times in 100 ul lower in 1.5 ml eppi
- dry in SpeedVac (35 C)

Chapter 13

Protocols - GC-ToF-MS measurement

13.1 Sample derivatisation

Materials:

- Methoxamine (MeOx)
- Pyridine (open under the hood only!)
- MSTFA
- Alkane mix (c10-c36) in Hexane
- chromacol vials and caps (big, small)

Mixtures:

- Solvent 1: 40 mg MeOx in 1 ml Pyridine
- Solvent 2: 10 ul Alkane mix in 1 ml MSTFA

Volumens of both solvents are shown for standard (small vol.) procedures.

Procedure:

- make sure samples are completely dry (1 h speed vac)
- add 20 ul (10 ul) of solvent 1 / sample
- incubate on rotary shaker, 30 C, for 60 min
- add 80 ul (25 ul) of solvent 2 / sample
- incubate on rotary shaker, 37 C, for 90 min
- centrifuge to spin down insoluble materials
- prepare aliquotes three times 28 ul or two times 15 ul (small glass vials)
- keep on room temperature until measurement (max. 10 days)

13.2 GC-MS measurement

needs to be written

Chapter 14

Data processing - MAUI

14.1 Processing In ChromaToF

Create a new folder in ChromaToF Pegasus Acquired Samples and import your files. The processing of files for Maui-assisted annotation is a two step process. Therefore two data processing methods have to be set up and applied to all files.

14.1.1 Resampling

Resampling is commonly applied and results into a data transformation enabling an improved detection of low abundant peaks and a reduction of noise. (Maybe include an example?)

The processing methods requires to tick **Export of** Subsequently, you are asked to define an output folder and the following paramter:

- Reduction rate: 4
- Beginning to end of the file
- .peg-files

14.1.2 Combo-export (.cdf & .csv)

Re-import the generated .peg-files into a subfolder and apply the following data processing method.

Activate the box **asddasd** and define for both file types the following parameter.

.cdf-file:

- export directory
-

.csv-file:

- export directory
-

14.2 Maui notes

14.3 Maui exports

With initiation of a project folder via `MTXQCvX2_init.Rmd` you created an `input`-folder containing three subfolders: `gc`, `inc`, `quant`. In the following all files that should be exported and copied into these folders are described in detail.

14.3.1 input/gc/...

Four input files are exported in order to assess the quality of the GC-MS performance of the run. The menu `Diagnostics` is selectable via right click on your Maui project name. Only the cinnamic acid peak areas are exported via the function `Export Quantification` with right click on the actual container.

Exported Files:

- `Alcane_intensities.csv` - Diagnostics/Export Alcane intensities
- `InternalStandard.csv` - cinAcid container, Export quantification, rename!
- `MassSum-73.csv` - Diagnostics/QC Mass Sum Export; enter: 73 for m/z 73
- `PeakDensities-Chroma.csv` - Diagnostics/ExportPeakDensities

14.3.2 input/quant/...

Only one container has to be exported and contains the peak areas of each metabolite and measurement. Keep in mind that you should have uploaded `pTop5 mass list` for the correct determination of peak areas in case of labeling with stable isotopes.

A further note - Maui performs absolute quantification and stores values in the `samplePeakGroups-QMQ` container. These quantities are determined by polynomial regression, and not linear regression.

The file `ManualQuantTable.tsv` is automatically generated by Maui during processing the absolute quantification.

Exported Files:

- `ManualQuantTable.tsv` - location: Maui-project/export/QM-AbsoluteQuantification/...¹
- `quantMassAreasMatrix.csv` - Quantification export of the container `samplePeakGroups`

14.3.3 input/inc/...

This exports are only required in case of an experiment including the application of stable isotopes. It's this the case you should have performed two things:

1. Used the optional upload of `pTop5 mass list`
2. Go through the pSIRM workflow in Maui

Exported Files:

- `DataMatrix.csv` - Export % Label of container `pSIRM-samplesPeakGroups`
- `pSIRM_SpectraData.csv` - pSIRM Spectra Export of container `pSIRM-samplesPeakGroups`²

¹Don't forget to rename it - e.g., e17123cz

²Requires the selection of `Natural_MIDs.txt`

Chapter 15

Data Processing - Metmax

15.1 Resampling

15.2 1D-basic

15.3 Reference search

15.4 Export for Metmax

15.5 Data extraction with Metmax

15.5.1 Peak areas

15.5.2 MIDs

Chapter 16

Frequently Asked Questions

16.1 What are additional quantification standards

16.2 How do I extend `conversion_metabolite.csv`

Chapter 17

Appendix - Tables

```
##
## Attaching package: 'dplyr'

## The following objects are masked from 'package:stats':
##
##   filter, lag

## The following objects are masked from 'package:base':
##
##   intersect, setdiff, setequal, union
```

This chapter shows the structure of all input or output csv-files that are referenced throughout the documentation. Please refer to the chapters for more detailed explanations.

17.1 MTXQC base tables

17.1.1 config_mtx tables

17.1.1.1 conv_filenames.csv

| Column.name | Description | Value |
|-------------------|---|-------------------------------------|
| Metabolite_manual | Manual defined metabolite name | #Alanine (2TMS) |
| Metabolite | Library name of the metabolite | Alanine_(2TMS)_BP_RI:1097_IDENT:B+C |
| Metabolite_short | Short version of library name of the metabolite | Alanine_(2TMS) |
| Lettercode | Lettercode version of metabolite name | Ala_2TMS |
| Q1_value | Checked if quant1:1 value available | x |
| Mass_Pos | m/z-value corresponding to m_inc | 118 |
| SE_sel | Evaluation of the MIDs | x |
| Q_sel | Evaluation for absolute quantification | x |
| nopsirm | Exclusively for absolute quantification | |
| Standards | Defined as standard | InternalStandard, Alk |

17.1.1.2 conv_filenames_manVal.csv

```
## 'data.frame':   3 obs. of  2 variables:
## $ AssociatedFile: Factor w/ 3 levels "inc","pSIRM_se",...: 3 2 1
## $ Filename      : Factor w/ 3 levels "DataMatrix_manVal.csv",...: 3 2 1
```

| | AssociatedFile | Filename |
|------|----------------|---------------------------------|
| 1 | sample_area | quantMassAreasMatrix_manVal.csv |
| 2 | pSIRM_se | pSIRM_SpectraData_manVal.csv |
| 3 | inc | DataMatrix_manVal.csv |
| NA | NA | NA |
| NA.1 | NA | NA |

17.1.1.3 DEPRICATED FILE: MQ_correction.csv

Not in use anymore since fluffy adventure!

```
## 'data.frame':   7 obs. of  3 variables:
## $ Metabolite : Factor w/ 7 levels "Alanine_(2TMS)_BP_RI:1097_IDENT:B+C",...: 1 2 3 6 7 4 5
## $ Cor_factor : num 12 12 6.33 11.82 2.41 ...
## $ Target_value: num 134695 134695 83259 388544 65154 ...
```

| Metabolite | Cor_factor | Target_value |
|-----------------------------|------------|--------------|
| Alanine_(2TMS)_BP_RI:1097_I | 12.00 | 134695.25 |
| Alanine_(3TMS)_MP_RI:1367_I | 12.00 | 134695.25 |
| Fructose_(1MEOX)(5TMS)_BP_ | 6.33 | 83259.33 |
| Glucose_(1MEOX)(5TMS)_MP_ | 11.82 | 388543.52 |
| Glycerol_(3TMS)_MP_RI:1280_ | 2.41 | 65153.65 |
| Fructose_(1MEOX)(5TMS)_MP_ | 6.33 | 83259.33 |
| Glucose_(1MEOX)(5TMS)_BP_] | 11.82 | 388543.52 |

17.1.2 config_files tables

These tables are supposed to be modified in relation to the individual needs of a project.

17.1.2.1 conversion_metabolite.csv

```
## 'data.frame':   166 obs. of  10 variables:
## $ Metabolite_manual: Factor w/ 115 levels "","#20G","#Adenine (2TMS)",...: 1 1 4 1 1 7 5 6 8 9 ...
## $ Metabolite       : Factor w/ 166 levels "Acetoacetic acid_(1MEOX)(1TMS)_MP_RI:1132_BMD",...: 1 2 4
## $ Metabolite_short : Factor w/ 161 levels "Acetoacetic acid_(1MEOX)(1TMS)",...: 1 2 4 5 7 6 8 9 10 1
## $ Lettercode       : Factor w/ 143 levels "2PGA","3PGA",...: 4 5 6 7 13 14 8 9 11 12 ...
## $ Q1_value         : Factor w/ 2 levels "","x": 1 1 1 1 2 2 2 2 1 1 ...
## $ Mass_Pos         : int NA NA NA NA NA NA NA 118 190 249 NA ...
## $ SE_sel           : Factor w/ 2 levels "","x": 1 1 1 1 1 1 2 2 2 2 ...
## $ Q_sel            : Factor w/ 2 levels "","x": 1 1 2 2 1 1 2 2 2 2 ...
## $ nopsirm          : Factor w/ 2 levels "","x": 2 2 1 2 2 2 1 1 1 2 ...
## $ Standards        : Factor w/ 3 levels "","Alk","InternalStandard": 1 1 1 1 1 1 1 1 1 1 ...
```

| Metabolite | Metabolite | Metabolite | Lettercode | Q1_value | Mass_Pos | SE_sel | Q_sel | nopsirm | Standards |
|---------------------|---------------------|-----------------|------------|----------|----------|--------|-------|---------|-----------|
| Acetoacet acid_(1M) | Acetoacet acid_(1M) | Aceto_1M | | | NA | | | x | |
| Acetoacet acid_(2T) | Acetoacet acid_(2T) | Aceto_2T | | | NA | | | x | |
| #Adenos(3TMS) | Adenosin_(3TMS) | Adenosin_(3TMS) | Adenosin_ | | NA | | x | | |
| | Adenosin_(3TMS) | Adenosin_(3TMS) | Adenosin_ | | NA | | x | x | |
| Alanine_beta_(2TMS) | Alanine_beta_(2TMS) | bAla_2T | x | | NA | | | x | |

17.1.2.2 letter_pathway_complete.csv

```
## 'data.frame': 100 obs. of 7 variables:
## $ Letter_Derivate: Factor w/ 69 levels "2HG","2OG","3PGA",...: 1 2 3 4 5 5 6 6 7 8 ...
## $ Lettercode : Factor w/ 100 levels "2HG","2OG","3PGA",...: 1 2 3 4 6 5 8 7 9 10 ...
## $ Pathway : Factor w/ 9 levels "aa","glut","glyc",...: 2 9 3 4 5 5 1 1 1 1 ...
## $ Pathway.1 : Factor w/ 9 levels "1-glyc","2-tca",...: 3 2 1 7 9 9 5 5 5 5 ...
## $ Met_pathway : Factor w/ 97 levels "1-glyc_DHAP_BP",...: 28 21 17 91 96 96 35 34 36 37 ...
## $ Subs_class : Factor w/ 11 levels "aa","amino sugar",...: 6 6 8 3 4 4 1 1 1 1 ...
## $ Met_class : Factor w/ 97 levels "aa_Ala_2TMS",...: 55 56 83 43 46 46 2 1 3 4 ...
```

| Letter_Derivate | Lettercode | Pathway | Pathway.1 | Met_pathway | Subs_class | Met_class |
|-----------------|------------|------------|--------------|----------------|--------------|----------------------------|
| 2HG | 2HG | glut | 3-glut | 3-glut_Glut_2h | organic acid | organic acid_Glut_2hydroxy |
| 2OG | 2OG | tca | 2-tca | 2-tca_Glut_2ox | organic acid | organic acid_Glut_2oxo |
| 3PGA | 3PGA | glyc | 1-glyc | 1-glyc_PGA | phosphate | phosphate_PGA |
| A | A | nucleobase | 7-nucleobase | 7-nucleobase_A | nucleobase | nucleobase_Adenosine |

17.1.2.3 incorp_calc_masses.csv

```
## 'data.frame': 90 obs. of 3 variables:
## $ Metabolite: Factor w/ 45 levels "Alanine_(2TMS)_BP_RI:1097_IDENT:B+C",...: 1 1 2 2 3 3 4 4 5 5 ...
## $ Mass_mz : int 116 118 188 190 245 249 273 275 273 276 ...
## $ LI_MID : Factor w/ 2 levels "m0","minc": 1 2 1 2 1 2 1 2 1 2 ...
```

| Metabolite | Mass_mz | LI_MID |
|-------------------------------------|---------|--------|
| Alanine_(2TMS)_BP_RI:1097_I | 116 | m0 |
| Alanine_(2TMS)_BP_RI:1097_I | 118 | minc |
| Alanine_(3TMS)_MP_RI:1367_I | 188 | m0 |
| Alanine_(3TMS)_MP_RI:1367_I | 190 | minc |
| Aspartic acid_(2TMS)_BP_RI:1433_IDE | 245 | m0 |

17.1.2.4 quant1_values.csv

```
## 'data.frame': 70 obs. of 3 variables:
## $ Letter_Derivate: Factor w/ 70 levels "2HG","2OG","3PGA",...: 1 2 3 4 5 6 7 8 9 10 ...
```

```
## $ Quant1_v4      : int  57270 34220 43480 7400 18710 134700 11480 22710 15030 11220 ...
## $ Quant1_v3      : int  57270 34220 43480 7400 18710 134700 11480 22710 15030 11220 ...
```

| Letter_Derivate | Quant1_v4 | Quant1_v3 |
|-----------------|-----------|-----------|
| 2HG | 57270 | 57270 |
| 2OG | 34220 | 34220 |
| 3PGA | 43480 | 43480 |
| A | 7400 | 7400 |
| Adenosine | 18710 | 18710 |

17.1.2.5 mid_backups.csv

```
## 'data.frame':    224 obs. of  4 variables:
## $ Metabolite      : Factor w/ 38 levels "Alanine_beta-(3TMS)_MP_RI:1435_IDENT:A+D",...: 1 1 1 1 2 2 ...
## $ Mass.m.z.       : int   188 189 190 191 116 117 118 119 188 189 ...
## $ BackupPeakArea: int   4960 876 307 53 2616179 323019 99834 19759 4960 876 ...
## $ BackupMID       : num   0.8005 0.1414 0.0495 0.0086 0.8553 ...
```

| Metabolite | Mass.m.z. | BackupPeakArea | BackupMID |
|--|-----------|----------------|-----------|
| Alanine_beta-(3TMS)_MP_RI:1435_IDENT:A+D | 188 | 4960 | 0.8005000 |
| Alanine_beta-(3TMS)_MP_RI:1435_IDENT:A+D | 189 | 876 | 0.1414000 |
| Alanine_beta-(3TMS)_MP_RI:1435_IDENT:A+D | 190 | 307 | 0.0495000 |
| Alanine_beta-(3TMS)_MP_RI:1435_IDENT:A+D | 191 | 53 | 0.0086000 |
| Alanine_(2TMS)_BP_R | 116 | 2616179 | 0.8552984 |
| Alanine_(2TMS)_BP_R | 117 | 323019 | 0.1056035 |
| Alanine_(2TMS)_BP_R | 118 | 99834 | 0.0326384 |
| Alanine_(2TMS)_BP_R | 119 | 19759 | 0.0064597 |
| Alanine_(3TMS)_MP_F | 188 | 4960 | 0.8005000 |
| Alanine_(3TMS)_MP_F | 189 | 876 | 0.1414000 |

17.2 Input data

17.2.1 MAUI derived tables

17.2.2 Metmax derived tables

17.3 Output data

17.3.1 Experimental Setup

17.3.2 MTXQCvX2 part1

17.3.2.1 output/gc/...

17.3.2.1.1 HM_GC_values.csv & qcmetric_xy.csv

MTXQC exports a file summarising quality factors for each of the four parameter evaluating the GC performance. A summary representing the values illustrated in the heatmap are shown in table `HM_GC_values.csv`, individual exports for each metric in table `qcmetric_xy.csv`.

```
## 'data.frame':   3 obs. of  3 variables:
## $ Column.name: Factor w/ 3 levels "Batch_Id","qc_metric",...: 1 2 3
## $ Description: Factor w/ 3 levels "Batch-Id","Class of QC metric",...: 1 3 2
## $ Value       : Factor w/ 3 levels "0.937254457",...: 3 1 2
```

| Column.name | Description | Value |
|-------------|--|-------------|
| Batch_Id | Batch-Id | e18274ba |
| qc_metric | QC metric factor corresponding with 1 - very good and 0 - very low | 0.937254457 |
| title | Class of QC metric | alkanes |

```
## 'data.frame':   3 obs. of  3 variables:
## $ Column.name: Factor w/ 3 levels "Batch_Id","qc_metric",...: 1 2 3
## $ Description: Factor w/ 3 levels "Batch-Id","Class of QC metric",...: 1 3 2
## $ Value       : Factor w/ 3 levels "0.937254457",...: 3 1 2
```

| Column.name | Description | Value |
|-------------|--|-------------|
| Batch_Id | Batch-Id | e18274ba |
| qc_metric | QC metric factor corresponding with 1 - very good and 0 - very low | 0.937254457 |
| title | Class of QC metric | alkanes |

17.3.2.1.2 IntStandard_normfactors.csv & IntStandard_stats.csv

```
## 'data.frame':   5 obs. of  3 variables:
## $ Column.name: Factor w/ 5 levels "Batch_Id","File",...: 2 5 1 4 3
## $ Description: Factor w/ 5 levels "Batch-Id","Determined normalisation factor",...: 4 5 1 2 3
## $ Value       : Factor w/ 5 levels "0.837457514",...: 4 2 3 1 5
```

| Column.name | Description | Value |
|-------------|---|-----------------|
| File | File name | e18274ba_17.cdf |
| PeakArea | Peak area of internal extraction | 89308492 |
| Batch_Id | Batch-Id | e18274ba |
| IntStd_fac | Determined normalisation factor | 0.837457514 |
| IntStd_eval | Evaluation of normalisation factor in relation to defined range plus/minus one standard deviation | within |

```
## 'data.frame':   8 obs. of  3 variables:
## $ Column.name: Factor w/ 8 levels "Batch_Id","File",...: 2 7 1 4 3 6 5 8
## $ Description: Factor w/ 8 levels "Batch-Id","Evaluation regarding QC",...: 3 7 1 5 2 6 4 8
## $ Value       : Factor w/ 8 levels "0.837457514",...: 7 5 6 1 8 4 2 3
```

| Column.name | Description | Value |
|-------------|---|-----------------|
| File | File name | e18274ba_17.cdf |
| PeakArea | Peak area of internal extraction standard | 89308492 |
| Batch_Id | Batch-Id | e18274ba |
| IntStd_fac | Normalisation factor | 0.837457514 |
| IntStd_eval | Evaluation regarding QC | within |

17.3.2.1.3 Min_Annotation.csv & SumArea_stats.csv

```
## 'data.frame': 9 obs. of 3 variables:
## $ Column.name: Factor w/ 9 levels "area_fac","Batch_Id",...: 3 2 6 9 7 4 8 1 5
## $ Description: Factor w/ 9 levels "Extracted Batch-Id derived from file name",...: 2 1 6 8 9 3 7 4 5
## $ Value : Factor w/ 9 levels "1.296568521",...: 9 8 2 6 3 5 4 1 7
```

| Column.name | Description | Value |
|-------------|---|-----------------|
| File | File name | e18274ba_17.cdf |
| Batch_Id | Extracted Batch-Id derived from file name | e18274ba |
| n_area | Number of peak areas per file | 101 |
| sum_area | Sum of all peak areas | 44614610885 |
| n_total | Total number of entries (including NA) | 107 |

```
## 'data.frame': 2 obs. of 3 variables:
## $ Column.name: Factor w/ 2 levels "Batch_Id","n_50": 1 2
## $ Description: Factor w/ 2 levels "Batch-Id","Number corresponding to fifty percent of the maximum number of annotated peaks per file": 1 2
## $ Value : Factor w/ 2 levels "53.5","e18274ba": 2 1
```

| | Column.name | Description | Value |
|------|-------------|---|----------|
| 1 | Batch_Id | Batch-Id | e18274ba |
| 2 | n_50 | Number corresponding to fifty percent of the maximum number of annotated peaks per file | 53.5 |
| NA | NA | NA | NA |
| NA.1 | NA | NA | NA |
| NA.2 | NA | NA | NA |

17.3.2.1.4 mz73_data.csv

```
## 'data.frame': 7 obs. of 3 variables:
## $ Column.name: Factor w/ 7 levels "Batch_Id","File",...: 2 1 3 6 4 7 5
## $ Description: Factor w/ 7 levels "Batch-Id","File name",...: 2 1 3 6 4 7 5
## $ Value : Factor w/ 7 levels "0.002407244",...: 7 6 3 2 4 5 1
```

| Column.name | Description | Value |
|-------------|--|-----------------|
| File | File name | e18274ba_17.cdf |
| Batch_Id | Batch-ID | e18274ba |
| mean_73 | Mean value of the sum of m/z 73 intensities per file | 16314646.1 |
| sd_73 | Standard deviation of the mean of the sum of m/z 73 intensities per file | 143890119.5 |
| n_peaks | Number of intensities used for statistics | 600 |

17.3.2.2 output/quant/...

17.3.2.2.1 calcheck_linearity.csv

```
## 'data.frame':   21 obs. of  3 variables:
##  $ Column.name: Factor w/ 21 levels "adj_r_squared",...: 9 2 4 17 14 21 10 11 7 16 ...
##  $ Description: Factor w/ 21 levels "", "Adjusted Rsquare value of linear regression of the calibrati
##  $ Value       : Factor w/ 18 levels "", "-898.3400476",...: 11 12 13 6 7 16 3 10 9 17 ...
```

| Column.name | Description | Value |
|-------------|--|-------------------------------------|
| Metabolite | Full library name of the metabolite | Alanine_(3TMS)_MP_RI:1367_IDENT:B+C |
| Batch_Id | Batch-ID | e18274ba |
| File | File name | e18274ba_53.cdf |
| QuantMasses | Defined quantification masses | 110.0 133.0 114.0 100.0 188.0 |
| PeakArea | Sum of peak areas based on defined QuantMasses | 190.0 12710956 |

17.3.2.2.2 CalculationFileData.csv

This is probably the most important file that is generated by running MTXQCvX2_part1.Rmd. It summarises all quality factors, experimental data and determined quantities of your experiment. This file provides the input for MTXQCvX2_part2-PostProcessing.Rmd.

```
## 'data.frame':   43 obs. of  4 variables:
##  $ Column.name: Factor w/ 43 levels "absconc","adj_r_squared",...: 4 6 9 13 36 41 42 12 43 23 ...
##  $ Class       : Factor w/ 9 levels "AnnExp","AnnExtract",...: 1 1 1 1 1 1 1 2 2 3 ...
##  $ Description: Factor w/ 41 levels "Absolute quantity in pmol",...: 3 14 14 16 6 14 40 7 9 18 ...
##  $ Value       : Factor w/ 37 levels "", "#Glycerol-3-phosphate (4TMS)",...: 27 24 25 28 1 3 34 18 26 2
```

| Column.name | Class | Description | Value |
|-------------|------------|---|-----------------|
| Batch_Id | AnnExp | Batch-Id extracted | e18274ba |
| CL | AnnExp | from file name Experimental | BE(2)-C |
| Cond | AnnExp | parameter Experimental | Control |
| File | AnnExp | File name | e18274ba_25.cdf |
| Standards | AnnExp | Defined as standard (InternalStandard, Alk) | |
| Time | AnnExp | Experimental parameter | 0 |
| Type | AnnExp | Type of measurement | sample |
| Extract_vol | AnnExtract | Defined extractes in count, mg or uL defined | 3290000 |
| Unit | AnnExtract | in Unit Defined unit for corresponding | count |
| Lettercode | AnnMet | Extract_vol Lettercode version of metabolite name | Glyc3P |

17.3.2.2.3 HeatMap_Quant_pTop5.csv

```
## 'data.frame': 5 obs. of 3 variables:
## $ Column.name: Factor w/ 5 levels "Batch_Id","Lettercode",...: 2 1 3 4 5
## $ Description: Factor w/ 5 levels "Batch-Id","Library name of metabolite",...: 4 1 2 3 5
## $ Value : Factor w/ 5 levels "0.996053496",...: 2 4 3 5 1
```

| Column.name | Description | Value |
|-------------|--|--|
| Lettercode | Short name of metabolite | Cit |
| Batch_Id | Batch-Id | e18274ba |
| Metabolite | Library name of metabolite | Citric acid_(4TMS)_MP_RI:1814_IDENT:B+D |
| Par | Parameter | R2_cal |
| Val | Value of the parameter for corresponding metabolite | 0.996053496 |

17.3.2.2.4 pTop5_Calibration_Samples_lincheck.csv

```
## 'data.frame': 7 obs. of 3 variables:
## $ Column.name: Factor w/ 7 levels "Batch_Id","count",...: 4 3 1 5 2 7 6
## $ Description: Factor w/ 7 levels "Batch-Id","Evaluation of peak area in relation to calibration cu
## $ Value : Factor w/ 5 levels "1","3PGA","51",...: 2 5 4 NA 3 3 1
```

| Column.name | Description | Value |
|-------------|--------------------------|-------|
| Lettercode | Short name of metabolite | 3PGA |

17.3.2.2.5 top5_CalibrationInfo_unique.csv

```
## 'data.frame': 8 obs. of 3 variables:
## $ Column.name: Factor w/ 8 levels "adj_r_squared",...: 6 5 2 7 1 4 8 3
## $ Description: Factor w/ 8 levels "Adjstuted Rsquare of calibration curve",...: 6 5 2 7 1 4 8 3
## $ Value : Factor w/ 8 levels "0.000194064",...: 6 5 7 8 2 4 1 3
```


| Column.name | Description | Value |
|---------------|---------------------------------------|---|
| Metabolite | Library name of metabolite | Citric acid_(4TMS)_MP_RI:1814_IDENT:B+D |
| Lettercode | Lettercode name of metabolite | Cit |
| Batch_Id | Batch-Id | e18274ba |
| Origin | Origin of quant1:1 value | Qstd |
| adj_r_squared | Adjusted Rsquare of calibration curve | 0.996053496 |

17.3.2.2.6 top5_QMQcurveInfo.csv

```
## 'data.frame': 15 obs. of 3 variables:
## $ Column.name: Factor w/ 15 levels "adj_r_squared",...: 8 7 14 10 2 5 3 4 13 11 ...
## $ Description: Factor w/ 15 levels "Adjusted Rsquare of calibration curve",...: 8 5 13 9 2 6 4 3 12 ...
## $ Value : Factor w/ 13 levels "0.000194064",...: 9 9 7 11 12 2 6 4 13 10 ...
```

| Column.name | Description | Value |
|------------------|---------------------------------------|---|
| Lettercode | Lettercode of metabolite name | Cit |
| Letter_Derivate | Derivate name | Cit |
| Quant1_v4 | Quant1:1 value in (pmol) | 52050 |
| Metabolite | Library name of metabolite | Citric acid_(4TMS)_MP_RI:1814_IDENT:B+D |
| Batch_Id | Batch-Id | e18274ba |
| Dilution | Dilution factor | 0.2 |
| ChromIntensities | Corresponding peak areas | 45074572 |
| Concentration | Concentration in (pmol) | 10410 |
| Origin | Origin of quantification standard | Qstd |
| Metabolite_short | Short name of metabolite | Citric acid 275_(4TMS) |
| adj_r_squared | Adjusted Rsquare of calibration curve | 0.996053496 |
| intercept | Intercept of calibration curve | 564.549288 |
| slope | Slope of calibration curve | 0.000194064 |
| max | Max. value of calibration curve | 52050 |
| min | Min. value of calibration curve | 260.25 |

17.3.2.3 output/inc/...

17.3.2.3.1 HeatMap_Incorporation.csv

```
## 'data.frame': 4 obs. of 3 variables:
## $ Column.name: Factor w/ 4 levels "Batch_Id","Lettercode",...: 2 1 3 4
## $ Description: Factor w/ 4 levels "Batch-Id","Lettercode name of metabolite",...: 2 1 3 4
## $ Value : Factor w/ 4 levels "0.740740741",...: 2 3 4 1
```

| | Column.name | Description | Value |
|----|-------------|---|-------------|
| 1 | Lettercode | Lettercode name of metabolite | 3PGA |
| 2 | Batch_Id | Batch-Id | e18274ba |
| 3 | Par | Parameter | NA_count |
| 4 | Val | Value of the parameter shown in heatmap | 0.740740741 |
| NA | NA | NA | NA |

17.3.2.3.2 SE_calculation_NAscore.csv

```
## 'data.frame':   5 obs. of  3 variables:
## $ Column.name: Factor w/ 5 levels "Batch_Id","fracr_prop",...: 3 1 5 4 2
## $ Description: Factor w/ 5 levels "Batch-Id","Class of NA-value",...: 4 1 2 5 3
## $ Value       : Factor w/ 5 levels "0","0.851851852",...: 4 5 1 3 2
```

| Column.name | Description | Value |
|-------------|---|-------------|
| Lettercode | Lettercode name of metabolite | Cit |
| Batch_Id | Batch-Id | e18274ba |
| na_frac_r | Class of NA-value | 0 |
| N | Number of MIDs | 46 |
| fracr_prop | Fraction of MIDs with x NA-values in relation to total number of evaluated MIDs | 0.851851852 |

17.3.2.3.3 SE_classification.csv

```
## 'data.frame':   9 obs. of  3 variables:
## $ Column.name: Factor w/ 9 levels "Batch_Id","count_score",...: 4 1 3 8 7 5 6 9 2
## $ Description: Factor w/ 6 levels "_expl_","Batch-Id",...: 4 2 3 1 1 1 1 6 5
## $ Value       : Factor w/ 9 levels "1861","2.30897367",...: 4 6 7 3 1 2 9 5 8
```

| Column.name | Description | Value |
|-------------|-------------------------------|-----------------|
| Lettercode | Lettercode name of metabolite | Cit |
| Batch_Id | Batch-Id | e18274ba |
| File | File name | e18274ba_26.cdf |
| sum_spa | _expl_ | 4297 |
| sum_bpa | _expl_ | 1861 |

17.3.2.3.4 SE_validation.csv

```
## 'data.frame':   6 obs. of  3 variables:
## $ Column.name: Factor w/ 6 levels "Batch_Id","count_score",...: 3 1 2 4 6 5
## $ Description: Factor w/ 6 levels "Batch-Id","Evaluation of MID",...: 3 1 2 5 4 6
## $ Value       : Factor w/ 5 levels "1","54","Cit",...: 3 4 5 2 2 1
```

| Column.name | Description | Value |
|-------------|-------------------------------|----------|
| Lettercode | Lettercode name of metabolite | Cit |
| Batch_Id | Batch-Id | e18274ba |
| count_score | Evaluation of MID | goodQ |
| N_count | Number of evaluated MIDs | 54 |
| sum_files | Number of detected MIDs | 54 |

17.3.2.4 MTXQCvX2 part3

17.3.2.5 MTXQCvX2 part4

Chapter 18

pSIRM experiments

The application of stable isotopes provides a powerful tool to track the activity of metabolic pathways. the time-dependent and atom-specific routing along a metabolic pathway resolved how substrates like glucose or glutamine are used in order to maintain a certain phenotype and energetic homeostasis.

We developed an approach called pulsed stable isotope resolved metabolomics (pSIRM) enabling the quantitative evaluation of metabolite pool sizes and incorporation of stable isotopes, e.g., $^{13}\text{C}_6$ -glucose. A thoughtful setup of the experimental design including the applied substrates and careful experimental handling are prerequisites for a successful pSIRM experiment. Essential aspects are collected in the below paragraphs along with a number of useful tweaks.

18.1 Experimental design

An *in vitro* pSIRM experiment lasts in total up to three days starting from the cell seeding at day zero. Further along the way up to two media changes should be included until the application of stable isotopes and harvesting the cells maintaining the continuous availability of nutrients and avoiding the accumulation of waste products (Figure ??fig:psirm)). The media change four hours prior the harvest is set up in order to give cells time to recover from the mechanical stress of the media change. At the time point of harvest cells should be in a perfect happy state regarding metabolic environment and stress.

Choose carefully the *seeding density of your cells* in the first place. High confluency inducing contact inhibition of cell growth has a strong impact on several cellular processes including the uptake of nutrients. Try to aim for petri dishes with a maximum confluency of 75-80 %. A pre-experiment including different cell densities for seeding at a number of experimental conditions helps you to get a feeling for the cell growth in general and an expected output of cells at the time point of the harvest. Later one is useful to plan sample extraction and measurement subsequently.

For adherent cell cultures only: Include for each experimental condition an additional petri dish that is solely used to determine the cell count at the time point of your harvest. This additional plate ensures a correct determination of absolute quantities and might reduce variation of pool sizes in the statistical analysis¹. Think carefully about control conditions and include cell culture dishes that are not labeled. These dishes function as a control for your labeling procedure and the natural abundance of isotopes.

¹Pelleting these cells and snap-freezing might give useful additional samples for western blotting.




Figure 18.1: Experimental design of a pSIRM experiment distinguishing short and long labeling with stable isotopes.

Table 18.1: Solvent composition of Full label medium (LM) and label buffer (BF) for a pSIRM experiment labeling with ^{13}C -glucose.

| Solvent | Base | Carbon source | Supplement |
|------------------------|--|--|-----------------|
| Full label medium (LM) | DMEM, without glucose, glutamine, pyruvate | ^{13}C -Glc (2.5 g/L) ^{12}C -Gln (2 mM) | small molecules |
| Label buffer (LB) | HEPES (5 mM), NaCl (140 mM), pH 7.4 | ^{13}C -Glc (2.5 g/L) ^{12}C -Gln (2 mM) | |

18.2 Experimental procedures

There is a slight difference in the protocols for adherent and suspension cells. Please read instructions and footnotes carefully.

Short² and long term labeling protocols differ only in the applied solvents during the labeling - either label buffer (LB) only or a combination of full label medium (LM) and label buffer (LB). Latter one is applied in order to remove extracellular metabolites of the media. Both label buffer and media contain the major nutrients / stable isotopes to keep the main substrates at constant supply at all times (Table 18.1). During the application of stable isotopes longer than a few minutes cells might sense the absence of further intermediates provided during standard cell culture procedure and adjust their metabolic program accordingly.

The quality of your data later heavily relies on the exact handling of the cells and a *consistent timing* throughout the pSIRM experiment. Especially the step removing the LB and quenching the cells should be a matter of a tenth of seconds rather than seconds. It is of great value to perform the cell harvest with a second person.

18.3 Protocol pSIRM

18.3.1 Adherent cell cultures

The herein described protocols are detailed explanations how to perform a pSIRM cell harvest for long term label application. If you want to label for less than 2 minutes omit solely *omit steps 6-8*.

Materials:

²I would rather recommend this up to 2 min of labeling

- Cell culture dishes, max. confluency 80 %
- Labeling media (LM) supplemented with substrates (5 ml / dish)³
- Label buffer (LB) supplemented with substrates (5 ml / dish)
- Ice-cold 50 % MeOH supplemented with 2 ug/ul cinnamic acid
- 2x 5 ml pipette and tips⁴
- Beaker
- Ice
- 15 ml falcons (chloroform resistant)
- Cell lifter
- Biological waste bin next to your bench

Procedure:

1. Pre-warm LB and LM in the water bath
2. Take a number of petri dishes (condition-wise including all biol. replicates)
3. Discard cell culture media
4. Carefully add *long term labeling* LM OR *short term labeling*: LB
5. Incubate cells on the bench or in an incubator
6. Discard LM (beaker)
7. Add immediatly 5 ml of LB
8. Rotate dish once in order to cover complete surface
9. Meanwhile 2nd person get prepared with 5 ml ice-cold MeOH
10. Discard LB into beaker and *immediatly* 2nd person quenches with ice-cold MeOH
11. Collect cell extracts using cell lifter
12. Transfer cell extracts into 15 ml falcons
13. Store falcons on ice until further processing

Repeat this procedure (step 6-10) for all dishes of a single condition first. Once MeOH is added metabolic processes are interrupted and cell extracts can be collected with the help of cell lifter without rush and subsequently transferred to 15 ml falcon and stored on ice until further processing (see chapter Cell extraction methanolic extracts).

Determine the cell count using your additional petri dishes for each condition.

18.3.2 Suspension cell cultures

Materials:

- Cell culture flasks
- Labeling media (LM) supplemented with substrates
- 5 ml pipette and tips⁵
- 1 ml pipette
- Beaker
- paper tissues
- Liquid nitrogen
- 15 ml falcons
- 1.5 ml eppendorf tubes
- Biological waste bin next to your bench

Procedure:

1. Pre-warm LM in the water bath
2. Determine the cell count of your cell suspension(s)

³Not required for short term labeling

⁴Highly recommended, makes labeling and harvest super quick

⁵Highly recommended, makes labeling and harvest super quick

3. Take aliquots of $10 - 15 \times 10^6$ cells and transfer into 15 ml falcon
4. Spin down cells very gently 300 g, 2 min at room temperature
5. Discard media into beaker
6. Resuspend cells gently in three-times 1 ml⁶
7. Incubate and keep warm
8. Fractionate cell suspension in three eppendorf tubes (3x 1 ml each)
9. Spin down quickly in top-bench centrifuge⁷
10. Discard media blandtly on paper tissues
11. Snap-freeze immediatly in liquid nitrogen
12. Store cells until further processing (see chapter Cell extraction suspension cells)

The important step here to be quick is the aliquotation of the cell suspension and subsequent spin down in the table centrifuge. Suspension cells are rather small, nevertheless 3×10^6 cells per extract are a good starting point for GC-MS measurements.

18.4 Hints & notes

- The only way to be reproducible and fast is to team up with a second person.
- Keep timing consistently through the experiment.
- Keep substrate concentrations constant throughout the experiment in all solutions.
- Supplement one stable isotopic labeled substrate with all remaining substrates in non-labeled form.
- Think about nutrient levels in your cell culture and your experimental conditions. Maybe you want to change things to physiological levels.
- Add additional plate to each condition in order to have material for western blotting and others.
- Check carefully the confluency of your dishes and determine seeding densities for different conditions.
- In case of small molecule inhibitors: Try to avoid to solve them in DMSO - strong impact on chromatography.

⁶To generate three replicates

⁷Most of the times 30 s are already enough

Chapter 19

GC-MS settings

In the following paragraphs details of GC-MS settings are described in detail. The herein described settings have been optimized for cell extracts measured in split-mode 1:5 on the instrument Pegasus 4D-C GC-ToF-MS in 1D mode equipped with an autosampler Gerstel MPS.

19.0.1 Autosampler settings

The table ?? summarises the defined parameter for method settings of the Gerstel MPS.

The injector itself is a temperature-regulated system providing the advantage of a focused injection of the sample. Gradient

19.0.2 Gas chromatography

The most important settings of the GC method are listed below. A graphical representation of the gradient in the oven is shown in figure 19.2 and corresponding values in table 19.1.

Flow path:

1. Inlet: Front
2. Capillary: GC Oven 50 m, 250 μ int. diameter, 0.25 μ film thickness, RTX-5 phase
3. Capillary: Detector 0.21 m, 250 μ int. diameter, 0.25 μ film thickness, RTX-5 phase
4. Detector: TOF

19.0.3 Mass spectrometer settings




Figure 19.1: Settings temperature-regulated injection for Gerstel MPS




Figure 19.2: GC gradient - graphical representation. Rate in ($^{\circ}\text{C}/\text{min}$), Target temperature in ($^{\circ}\text{C}$), Duration in (min).

Table 19.1: GC gradient profile - Rate shown in (K/min), Target Temp in ($^{\circ}\text{C}$), Duration in (min)

| Rate | Target_Temp | Duration |
|---------|-------------|----------|
| Initial | 68 | 2 |
| 5 | 120 | 0 |
| 7 | 200 | 0 |
| 12 | 320 | 6 |

Chapter 20

Solutions and Materials

20.1 Solutions

20.1.1 Label buffer

Materials:

- ddH₂O (500 ml)
- 140 mM NaCl (4.1 g)
- 5 mM Hepes (0.569 g)
- pH calibration 7.4

Procedure:

- Weigh the correct amounts of Hepes and NaCl
- Resolve in a glass bottle with 450 ml of water
- Stir carefully
- Check and adjust pH
- Adjust volumes to 500 ml

20.1.2 MCW

Materials:

- Methanol
- Chloroform
- ddH₂O
- Cinnamic acid stock in MeOH (2 mg/ml): final conc. 2 µg/ml

Procedure:

- Mix the solvents in the ratio of volumes - Methanol:Chloroform:Water – 5:2:1
- Supplement cinnamic acid stock 1:1000
- Store at -25°C

20.1.3 Alkane-Mix

Materials:

Table 20.1: List of chemicals including company names and catalog number.

| Chemical | Catalog.Nb | Company | Storage |
|--|--------------|----------------|---------|
| Chloroform | 132950-1L | Sigma | RT |
| Cinnamic acid | C80857 5g | Aldrich | RT |
| Hexane for HPLC > 97% (GC) | 34859-1L | Sigma | RT |
| Methanol | 1060351000 | Merck | RT |
| Methoxyamine hydrochloride | 226904-5G | Aldrich | RT |
| MSTFA - Methyltrimethylsilyltrifluoroacetamide | 701270201 | Macherey-Nagel | 4°C |
| Pyridine | 270970-100ML | Sigma | RT |

- Hexane
- Alkanes: c10, c12, c15, c17, c19, c22, c28, c32, c36
- Thermo mixer
- Glass vials and caps

Procedure:

- Prepare stock solutions in hexane:
- c10 - c17 (liquid): 25 ul/ml
- c19 - c32: 20 mg/ml
- c36: two-times 15 mg/1.5 ml
- Warm up alkane stocks in thermo mixer 40°C
- Prepare a test mixture in equal amounts, e.g., 50 ul each, but use twice the volume of c36
- Mix test mixture with MSTFA: 10 ul / 1 ml MSTFA
- Check alkane profile by GC-MS
- If required: adjust volumes and re-test or create larger volume of your mixture for aliquots
- Store aliquots in glass vials, close well and store at 4°C
- For usage: gently warm up glass vials at 30 C on thermo mixer for 10 min and vortex before adding it to the MSTFA

Adjust the volumes of the alkane stocks in order to create a curve shaped distribution of all alkanes in the chromatogram: lower intensities for c10 and c32-36, slowly increasing intensities for the alkanes in between.

20.2 Materials

Bibliography

Pietzke, M., Zasada, C., Mudrich, S., and Kempa, S. (2014). Decoding the dynamics of cellular metabolism and the action of 3-bromopyruvate and 2-deoxyglucose using pulsed stable isotope-resolved metabolomics. *Cancer & metabolism*, 2(1):9.